

2. pSC101-dCas9-eGFP-M library

2.1 Obtain pSC101-dCas9-eGFP-M-1 fragment by PCR

The target fragment is located on the pSC101-dCas9-eGFP plasmid, and the amplified target fragment is obtained by PCR. The PCR system and procedure are showed below.

PCR system (50 μ L)	
pSC101-dCas9-eGFP	10 ng
M1-F=0719	2 μ L
1-R	2 μ L
2 x Mix	25 μ L
DDW	To 50 μ L

Table 1

The PCR products are detected by agarose gel electrophoresis, and the correct target fragment is 4393bp. We obtained the correct target fragment, and the sample is purified.

2.2 Obtain pSC101-dCas9-eGFP-M-2 fragment by PCR

The target fragment is located on the pSC101-dCas9-eGFP plasmid, and the amplified target fragment is obtained by PCR. The PCR system and procedure are showed below.

PCR system (50 μ L)	
pSC101-dCas9-eGFP	10 ng
2-F	2 μ L
M2-R	2 μ L
2 x Mix	25 μ L
DDW	To 50 μ L

Table 2

The PCR products are detected by agarose gel electrophoresis, and the correct target fragment is 4820bp. We obtained the correct target fragment, and the sample is purified.

2.3 Gibson connection

The pSC101-dCas9-eGFP-M-1 fragment and pSC101-dCas9-eGFP-M-2 fragment are connected by Gibson connection method, and the connection system is as follows.

Connection system (10 μ L)	
pSC101-dCas9-eGFP-M-1	1 μ L
pSC101-dCas9-eGFP-M-2	1 μ L
2 x ClonExpressMix	5 μ L
DDW	3 μ L

Table 3

2.4 Colony PCR

After the petri dish is incubated at 37°C for 12 hours, 30 colonies were selected on the plate. The colony PCR system and procedure were as follows.

PCR system (10 μL)	
pSC101-dCas9-eGFP-M	1 μ L
link-GG-F=0615	0.2 μ L
link-GG-F=0615	0.2 μ L
2 x Mix	5 μ L
DDW	3.6 μ L

Table 4

The correct target fragment is about 1711bp, and the positive rate was approximately 60%. We collectively transferred the positive monoclonal colonies identified through colony PCR into a single LB medium for Extended cultivation, and the plasmids were put forward as a mutation library.